



# Clearing the TEMP Folder

## TMP and TEMP



- Close all applications and make sure your files are saved
- Open **File Explorer**
- Enter %TMP or %TEMP% in the **Address Bar**
- This should take you to a file path similar to C:\Users\Scott\AppData\Local\Temp
- Remove the data in this folder (or rename it to Temp2 if you want to keep it around)

## Windows Temp



- Go to C:\Windows\Temp to remove the global temp folder, if that is needed

## Caution: Open Apps

**Note:** deleting temp files is safe, but if Windows OS shows a popup indicating that the file is in use, then do not delete that file because an open app is using that temp file

From:

<https://help.go-ci.com/> - **Crystal Instruments Help**

Permanent link:

<https://help.go-ci.com/general:temp-folder?rev=1718732363>

Last update: **2024/06/18 17:39**